

Product tips

Aroma base and Aroma sticks

Requirements:

Aroma base	ZAT-02001
Perfume oil	YZS-xxxxx
Measuring beaker	
Spoon	
Glass bottle, e.g.	ZVE-00011
Shrink capsule, e.g.	ZVE-70011
Hot air gun	FHN-00001
Aroma sticks	ZAT-01001 (Ø 2,5 mm) or ZAT-01002 (Ø 3 mm)



Recipe for 100 ml:

Aroma base 90 ml

Perfume xxx 10 ml

Mix the two ingredients thoroughly and immediately pour the mixture into a glass bottle. Find a suitable bottle in the Gildewerk web shop at www.shop.gildewerk.com.

For extra security, the bottle and cork can be sealed with an appropriate shrink capsule. Cut the capsule to the correct length - the cork and part of the bottle's neck must be covered. By briefly heating this with the hot air gun, the capsule will shrink tightly around the cork/bottle.

Cut the sticks to the desired length (max. 35 cm). These are available in two thicknesses, and whichever you use is your choice.

For effective fragrance distribution, we recommend that you use at least 5 sticks. Our special Aroma base guarantees even and long-lasting fragrance distribution.

Aroma base and Aroma sticks

- Aroma sticks are great for use in bedrooms, bathrooms and guarantee a long lasting fragrance for the whole living area.
- Arrange the natural sticks in the fragrance filled bottle.
- To optimize the aroma, flip the sticks over once a week. (wash hands afterwards with water and soap)
- Remove spilled material immediately to avoid spots and discolouration. **Warning:**
- The material is flammable, do not use presence of open fire, don't smoke.



- *use this symbol on your label*
- Keep out of reach of children.